

# A site for sore eyes.



100% recycled and recyclable.

Easier access.

Larger, shatterproof mirror.

50% smaller - less plastic.†

Illustrated instructions, so no translation required.

Mirror is 88% bigger.†

Uses 47% less wall space.†

## Reliwash™

### EYE WASH STATION

Now made with **100% recycled and recyclable plastic** this updated station is manufactured in high visibility green from recycled plastic materials incorporating two 500ml bottles of sterile solution, two eye pads, and an extra-large shatterproof mirror. With isolated workers in mind, easy to follow instructions on the fascia make self-administration by an injured person practical and straightforward.

The innovative wall mounting layout is almost half the size of traditional eye wash stations, and now makes it possible to deliver a complete response for emergency eye injuries in virtually any location.

Size with bottles: 19cmW x 30.7cmH x 5.9cmD

Size without bottles: 19cmW x 15.8cmH x 5.9cmD

906: EYE POD STATION - STANDARD BOTTLE  
CONTENTS: 2 x Dressings, Eye Pad, 2 x Eye Wash, Bottles 500ml, 1 x Mirror

909: EYE POD STATION - EYE CAP BOTTLE  
CONTENTS: 2 x Dressings, Eye Pad, 2 x Eye Wash, Bottles 500ml, 1 x Mirror

\* Excludes label and product packaging.  
† Compared to the previous design.

**Better for the workplace, better for the planet.**