



aura³

Aura 3 is the latest version of our industry leading first aid kit box.



KEY UNIQUE FEATURES:

New modern design.

An industry first: 100% recyclable and made from 100% recycled material.

No Labels Boxes are embossed and screen printed on the lid, laser-etched on the base.

Built-in mounting bracket

Easy one-handed opening when first aid is needed in a hurry.

WHY IS AURA3 BETTER FOR THE WORKPLACE AND THE PLANET?

Aura3 represents a first aid industry first. Aura3 is the first and only first aid kit box made from 100% recycled plastic and 100% recyclable.

Based on weight, Aura3 will save approx. 3.42 million plastic bottles each year. This many plastic bottles would cover the pitch at Wembley twice.

DOES AURA3 USE LABELS?

Labels are no longer used on Aura. Each Aura3 box has **embossed and screen-printed graphics** on the lid to denote its use – e.g. First Aid Kits carry an embossed cross, Eyewash kits feature an embossed eyewash symbol etc.

Multilingual design. By using symbols, Aura3 kits are easily recognisable, and each kit's intended usage can be easily understood in all languages. Other graphics are screen printed on the lid. E.g. Workplace kits carry the reliance logo, burns kits carry the reliance logo and a burns icon. Boxes are **laser printed on the base** – this reduces the use of plastics even more as labels are not required.

WHAT ELSE IS DIFFERENT ABOUT AURA 3?

Aura3 features a single clip that securely holds the lid in place – this means the box can quickly and easily opened with one hand. The box can also be secured with an anti-tamper tag to ensure contents are always complete.

CAN AURA3 BE MOUNTED TO WALL?

Aura features an integral fixing bracket that can be screwed to the wall - allowing the Aura box to slide in place over the bracket and be securely fastened to the wall. The bracket in Aura3 is smaller and neater than previous mounting brackets.



Better for the workplace, better for the planet.





WHAT SIZES ARE AVAILABLE IN AURA3?

Aura is available in large, medium, small and extra small sizes. The previous compact aura range no longer exists - Aura3 extra small is the previous large compact Aura in an updated design that complements the Aura3 range – featuring the same embossed graphics and laser etched base.

Aura3 large - 30.5cmH x 34.5cmW x 10.5cmD

Aura3 medium – 25.5cmH x 28cmW x 10cmD

Aura3 small – 21cmH x 23.5cmW x 10cmD

Aura3 extra small – 20.5cmH x 20cmW x 6.5D

What Colours are available?

Aura3 boxes feature glossy, eye catching colour coded lids – available in Green, Blue, Yellow, Red and Orange (see our range leaflet for full details).

ARE AURA3 KITS COMPLIANT?

Aura3 Workplace kits meet the BS8599-1:2019 standard HSE versions of the kits are also available.

BS8599-1:2019 first aid kits meet a higher standard than HSE kit requirements including the right products to comply with health and safety regulations. The new standards (BS 8599-1:2019) came into effect on 10th January 2019. The new Workplace First Aid Kits contain more disposable gloves, gel burns dressings, fewer medium sterile dressings, but more large sterile dressings, foil blankets, safety scissors and fewer triangular bandages.

ARE AURA3 EMPTY BOXES AVAILABLE?

Whilst we are aiming to only bring in complete kits from RMSH, a supply of Aura3 empty boxes will be available in small, medium and large sizes in green and extra small in orange. Empty boxes are multi use and feature the first aid cross only on the lid.

CAN CUSTOMERS HAVE THEIR OWN BRAND VERSION OF AURA3?

Customers can have the option to have their own Aura3 range of kits – by having their logo screen printed on the lid.

This option is available on a minimum order quantity of 1000 boxes.

Better for the workplace, better for the planet.