

\* Excludes product and packaging.



Integral bracket





**Embossed cross** 



Easy identification



Tamper seal



Laser etching



# the right kit for your workplace.

### **BS8599-1:2019 COMPLIANT WORKPLACE FIRST AID KITS**

Our New Aura 3 BS 8599-1:2019 kits are better for the workplace and better for the planet, providing quick and easy access to first aid, fully complaint with the latest standard requirements. Aura3 represents a first for the first aid industry, the first 100% recycled and fully recyclable first aid kit box on the market.

The fresh new box design feature colour coded, embossed, and screen-printed lids. The new graphics stand out from a distance, enabling quick identification in an emergency, especially in a busy workplace. Compact and convenient, the integral bracket allows easy placement and fixing.

British Standard compliant first aid kits reflect the most up-to-date changes in both workplace practice and risks, making them demonstrably more 'fit for purpose' in today's environment.

The quantities and sizes of first aid kits to be used in a workplace environment should be determined by the number of people on site and the category of hazard they are exposed to.

Better for the workplace, better for the planet.



# What does the law say?

The Health and Safety (First Aid) regulations 1981 states:

"An employer shall provide or ensure that there are provided such equipment and facilities as are adequate and appropriate in the circumstances for enabling first aid to be rendered to his employees if they are injured or become ill at work"

The following information categorises workplaces by level of risk, and then recommends an appropriate size of kit, depending on the number of employees and members of the public exposed to those risks.

#### DO A RISK ASSESSMENT.

LOW RISK: (e.g. shops, offices, libraries etc.) HIGH RISK: (e.g. light engineering and assembly work, food processing, warehousing, etc.)

# DECIDE ON THE SIZE AND QUANTITY OF KITS REQUIRED TO MEET YOUR RISK NEEDS ASSESSMENT.

SMALL KIT USAGE GUIDELINES

LOW RISK: 1 kit for less than 25 employees HIGH RISK: 1 kit for less than 5 employees

MEDIUM KIT USAGE GUIDELINES

LOW RISK: 1 kit per 25-100 employees HIGH RISK: 1 kit per 5-25 employees

LARGE KIT USAGE GUIDELINES

1 kit per 100 employees LOW RISK: HIGH RISK: 1 kit per 25 employees

# TRAVEL AND MOTORING

Staff working remotely from a company vehicle should be equipped with a travel and motoring first aid kit. This does not supersede BS8599-2 for vehicle first aid kits, employers should refer to BS8599-2 when equipping vehicles with first aid kits.

# 1 Medium Kit

	330 SMALL	343 MEDIUM	348 LARGE	3011 TRAVEL &
				MOTORING <sup>1</sup>
Bandage, Conforming	1	2	2	0
Bandage, Triangular	2	3	4	1
Dressing Pad, Adhesive	0	0	0	1
Dressing, Burn <sup>3</sup>	1	2	2	2
Dressing, Eye Pad	2	3	4	0
Dressing, Finger	2	3	4	0
Dressing, Sterile, Large	2	3	4	0
Dressing, Sterile, Medium	2	4	6	1
Dressing, Trauma, Medium	0	0	0	1
Foil Blanket	1	2	3	1
Gloves, Nitrile (Pairs)	6	9	12	2
Guidance Leaflet	1	1	1	1
Plasters, Washproof, Assorted	40	60	100	10
Resuscitation Face Shield <sup>2</sup>	1	1	2	1
Shears	1	1	1	1
Tape, Microporous	1	2	3	0
Wipes	20	30	40	10

<sup>1</sup> The quantities in the travel and motoring first aid kit and personal issue first aid kit are appropriate for one person only. 2 The shield is used to protect a first-aider from infection from a casualty during mouth-to-mouth resuscitation. 3 A Burn Dressing cools the burn injury without any prior treatment with water or any other medium. The dressing does not dry out within one hour of application.



#### **BS8599-1:2019 WORKPLACE**



SIZE: UNIT:

MOQ:

343

SIZE:

UNIT: MOQ:

348

SIZE:

MOO:

3011

BS8599-1:2019 SMALL WORKPLACE KIT in Small Green Aura3 Box 21cmH x 23.5cmW x 10cmD 10



BS8599-1:2019 MEDIUM WORKPLACE KIT in Medium Green Aura3 Box 25.5cmH x 28cmW x 10cmD 10



BS8599-1:2019 LARGE WORKPLACE KIT in Large Green Aura3 Box 30.5cmH x 34.5cmW x 10.5cmD

#### BS8599-1:2019 & BS8599-2:2014 MEDIUM TRAVEL & MOTORING KIT



BS8599-1:2019 & BS8599-2:2014 COMPLIANT TRAVEL & MOTORING KIT in Xtra Small Green Aura3 Box

SIZE: 20.5cmH x 20cmW x 6.5cmD 20

The BS8599-1:2019 Guidelines recommend supplementary kits as an extra resource to British Standard Compliant Small, Medium, Large and Travel & Motoring First Aid Kits.

### PERSONAL ISSUE FIRST AID KIT

Containing basic first aid components for one person. e.g. maintenance, security personnel, cleaners, first aiders, etc.

#### **CRITICAL INJURY PACK**

Designed for environments where people could sustain a severe injury. Contains items to help manage catastrophic blood loss. e.g. work with dangerous machinery or sharp instruments, cutting equipment, power tools, construction, agriculture, forestry, etc.



