

Hold back the burn

NEW
INNOVATIVE
Glove Design

CoolTherm™ Professional Emergency Dressing for Burns

ISOLATES each finger to reduce recovery time and damage to skin - cooling the heat and easing the pain



CoolTherm™
cools the heat - eases the pain
Professional Emergency Dressing for Burns



+44 (0) 08456 448808 | sales@reliancemedical.co.uk | www.reliancemedical.co.uk

RELIANCE MEDICAL
'Leading First Aid'



REV61223

CoolTherm™

cools the heat - eases the pain

Professional Emergency Dressing for Burns

Keep on hand to treat burns

- Sterile, non-woven trauma hydrogel burn dressing
- Minimises trauma and skin damage
- Non-adherent to burnt skin
- 5 years shelf life



COOLTHERM HAND DRESSING

5926 COOLTHERM HAND DRESSING

SIZE: Hand
UNIT: 1
MOQ: 10



Burns to hands present additional complications. The skin between fingers can adhere together, leading to webbing. This can result in a long and complex recovery, including skin grafts. Avoiding this is crucially important in the first aid treatment of burns to hands.

The CoolTherm™ 'Glove' hand dressing (**fig.1**) is a unique design, made from a non-woven carrier dressing that is very soft, gentle and highly conformable.

The 'Glove' includes exclusive finger separators to prevent webbing (**fig.2**).


The gel saturated mitt is easily donned, promptly covering all digits with soothing, cooling gel, and providing immediate analgesic effects to the casualty.

The fingertips are easily accessible for observation (**fig.3**), as are the thumb and small finger. Allows pinching of fingernail to check healthy blood flow. The mitt can be easily removed, by lifting it back, or cutting away.




HEAT EXCHANGE

CoolTherm™ absorbs the thermal energy from the burn, dissipating the heat, and cooling the skin...



PAIN RELIEF

CoolTherm™ relieves the pain by isolating damaged nerve endings from the air...



BARRIER TO CONTAMINANTS

CoolTherm™ protects against contamination by providing a sterile physical barrier...

CoolTherm™
cools the heat - eases the pain
Professional Emergency Dressing for Burns

A New Generation
of Professional
Emergency Burn
Dressings have arrived.