



**hypothermia blankets**

COVER ↻ PROTECT ↻ SURVIVE








The **Igloo Hypothermia Blanket** helps to preserve lives during initial medical care. The Igloo blanket combines state-of-the-art materials and technology to **regulate temperature swiftly and effectively** during emergency situations, playing a critical role in initial medical assessment and interventions.

Designed for emergency situations, The Igloo Blanket combines **cutting-edge technology** with practical features to quickly stabilise temperature and maintain optimal thermal retention.

The Igloo Blanket is the ultimate choice for the Military and Emergency Services including dedicated teams trained to save lives in rugged terrain. Available in **Emergency Orange** and **Coyote Brown** for military use.





# Life-saving heat to the core

Why is the Igloo Blanket a survival essential?




## Effective thermal stabilisation

**Ultra-light, ultra effective:** The Igloo Blanket excels in thermal retention and heat distribution, even under extreme conditions. This remarkable insulation capability is made possible by Igloo technology, an ultra-light fabric known for its exceptional ability to retain heat even under extreme conditions.

**Minimising heat loss:** The blanket features an innovative network of cells that capture warm air - allowing it to flow.

**Reflecting the heat back:** A reflective metallic layer directs heat back toward the user, enhancing insulation and locking in body heat.



**LIFE-SAVING HEAT:**  
CELLS TRAPPING AIR

**ELASTIC:**  
HUGS BODY

**REFLECTIVE COATING:**  
REFLECTS HEAT

## 2 Rapid deployment

In emergencies, every second counts. The Igloo Blanket unfolds quickly and effortlessly and starts to maintain body temperature in seconds. Due to its small space-saving pack, the Igloo Blanket is an essential for all first responders and disaster relief teams.

Igloo Blankets are not pressure packed and can be opened, inflated and deployed more rapidly as a result.



**EASILY ADAPTS  
TO USER'S SHAPE**



COVER  
PROTECT  
SURVIVE





**Body-hugging fit:** The Igloo Blanket adapts to the user's shape, snugly wrapping them. Its cocooning effect means no gaps and no heat escape, but is adaptable if a looser fit is required.

**Safe for all:** The Igloo Blanket is latex free, ensuring compatibility with all patients and compliance with hospital entrance safety protocols.

**Optimised heat distribution:** The Igloo Blanket quickly circulates and maintains warmth evenly across the body.



**Easy to seal and reopen.**



### Reliable thermal performance from incident to ER

**Trauma management:** Right from the scene of the incident, from prehospital to the operating theatre, the Igloo blanket maintains seamless life-saving heat - maximising survival outcomes.

**Radiation and conduction minimisation:** The technology reflects body heat back to the user, maintaining core temperature.

**Locks in heat.** From climbing accident to military medical extraction, the impermeable design shields against wind-chill, blocks evaporation and prevents heat loss whatever the environment.



## It's a wrap...

Compact, ultra-light, ultra effective.

**Igloo is the final word in life-saving thermal retention.**



# What blanket should I use?



## 2 layer hypothermia blanket

A 2-layer blanket:  
Designed to maintain body heat.



**IG010**  
IGLOO CORE 2-LAYER  
HYPOTHERMIA BLANKET  
**MEDIUM | Emergency Orange**  
2.3m x 1.5m | 5 TOG



**IG009**  
IGLOO CORE 2-LAYER  
HYPOTHERMIA BLANKET  
**MEDIUM | Coyote Brown**  
2.3m x 1.5m | 5 TOG



**IG012**  
IGLOO PRO  
HYPOTHERMIA BLANKET  
**SMALL | Emergency Orange**  
1.1m x 1.5m | 5 TOG



**IG011**  
IGLOO CORE 2-LAYER  
HYPOTHERMIA BLANKET  
**SMALL | Coyote Brown**  
1.1m x 1.5m | 5 TOG





# PRO

## 3 layer hypothermia blanket

A 3-layer blanket: Designed to maintain body heat in scenarios involving critical, life-threatening, injuries and conditions.



**IG004**  
IGLOO PRO 3-LAYER HYPOTHERMIA BLANKET  
**X LARGE | Emergency Orange**  
2.3m x 2m | 7 TOG



**IG003**  
IGLOO PRO 3-LAYER HYPOTHERMIA BLANKET  
**X LARGE | Coyote Brown**  
2.3m x 2m | 7 TOG



**IG002**  
IGLOO PRO 3-LAYER HYPOTHERMIA BLANKET  
**LARGE | Emergency Orange**  
2.3m x 1.8m | 7 TOG



**IG001**  
IGLOO PRO 3-LAYER HYPOTHERMIA BLANKET  
**LARGE | Coyote Brown**  
2.3m x 1.8m | 7 TOG



**IG008**  
IGLOO PRO 3-LAYER HYPOTHERMIA BLANKET  
**MEDIUM | Emergency Orange**  
2.3m x 1.5m | 7 TOG



**IG007**  
IGLOO PRO 3-LAYER HYPOTHERMIA BLANKET  
**MEDIUM | Coyote Brown**  
2.3m x 1.5m | 7 TOG



**IG015**  
IGLOO PRO HYPOTHERMIA BLANKET  
**BABY | Emergency Orange**  
1.1m x 1.5m | 7 TOG

# XTREME

## 3-layer hypothermia blanket + heat

A 3-layer blanket with heat pads: Designed to maintain body heat in scenarios involving critical, life threatening, injuries and conditions.



**IG006**  
 IGLOO XTREME 3-LAYER HYPOTHERMIA BLANKET + HEAT  
 X LARGE | Emergency Orange  
 2.3m x 2m | 7 TOG



**IG005**  
 IGLOO XTREME 3-LAYER HYPOTHERMIA BLANKET + HEAT  
 X LARGE | Coyote Brown  
 2.3m x 2m | 7 TOG

with  
**active heat**







# It's a wrap...

Compact, ultra-light, ultra effective.

**Igloo is the final word in life-saving thermal retention.**



Unit 26  
Llandygai Industrial Estate  
Llandygai  
Bangor  
Gwynedd  
LL57 4YH  
T: +44 (0)8456 448808

CONTENT SUBJECT TO CHANGE.